











Screw Cap Vials for Compatible Systems

	Clear	Amber	Amber Max Recovery	Clear Glass Max Recovery	Qsert Vial	Amber Qsert	PP 300 µL	PP 750 µL	10 mm Cap	PP 250 µL 8 mm Cap
12 × 32 mm										
Vial Number	50	51	52	53	54	55	56	57	58	59
Compatible Systems										
Agilent Technologies	▪	▪	▪	▪	▪	▪	▪	▪	—	—
Alcott, Antek, CTC, Spark Thermal Separations	—	—	—	—	—	—	—	—	▪	▪
Beckman, Dynatech, Finnigan, Fisons, Gilson	▪	▪	▪	▪	▪	▪	▪	▪	—	—
Hitachi, LDC, Perkin-Elmer, Shimadzu, Spectra-Physics, Thermo, Varian	▪	▪	▪	▪	▪	▪	▪	▪	▪	▪
Cap and Septum										
PE Septumless Caps	186004169	186004169	186004169	186004169	186004169	186004169	186004169	186004169	—	—
Black Cap	—	—	—	—	—	—	—	—	WAT058875	186004717
Cap and Septum, Silicone/ PTFE, Assembled	—	—	—	—	—	—	—	—	—	WAT094174
Septum Only, PTFE/ Silicone, Pre-slit	—	—	—	—	—	—	—	—	—	WAT058876
Septum Only, Silicone/PTFE	—	—	—	—	—	—	—	—	WAT058874	WAT210685
Septum Only, PTFE	—	—	—	—	—	—	—	—	—	WAT058886
Screw Cap and Septum-Silicone/PTFE										
Blue LectraBond	186000274	186000274	186000274	186000274	186000274	186000274	186000274	186000274	—	—
Red LectraBond	186002129	186002129	186002129	186002129	186002129	186002129	186002129	186002129	—	—
Green LectraBond	186002130	186002130	186002130	186002130	186002130	186002130	186002130	186002130	—	—
Screw Cap and Pre-slit Septum-Silicone/PTFE										
Blue LectraBond	186000305	186000305	186000305	186000305	186000305	186000305	186000305	186000305	—	—
Red LectraBond	186002128	186002128	186002128	186002128	186002128	186002128	186002128	186002128	—	—
Green LectraBond	186002127	186002127	186002127	186002127	186002127	186002127	186002127	186002127	—	—
Storage Cap										
Black Solid Cap 9 mm with Silicone/PTFE Liner	186007187	186007187	186007187	186007187	186007187	186007187	186007187	186007187	—	—

All items come in quantities of 100 unless otherwise noted.