




## DisQuE EXTRACTION AND CLEANUP TUBES AND POUCHES

DisQuE Extraction and Cleanup Tubes and Pouches are available separately for customized applications and method development. The salts contained in the 50 mL tubes are also available in a pouch format for greater flexibility. The cleanup tubes are available in a standard 2 mL size as well as a 15 mL size for sample enrichment.

### Ordering Information

#### DisQuE Dispersive Sample Preparation Products

Description	P/N
<b>Individual Extraction Tubes (Tube 1)</b>	
50 mL Empty Tube for QuEChERS Extraction	186006814
DisQuE 50 mL Tube/ AOAC-Acetate	
 DisQuE 50 mL tube containing: 1.5 g Sodium Acetate and 6 g MgSO <sub>4</sub>	186004571
DisQuE 50 mL Tube/ CEN-Citrate	
 DisQuE 50 mL tube containing: 1 g trisodium citrate dihydrate, 0.5 g disodium hydrogencitrate sesquihydrate, 1 g NaCl and 4 g MgSO <sub>4</sub>	186004837

Description	P/N
<b>Individual Extraction Pouch</b>	
DisQuE Pouch	
 1.5 g sodium acetate, 6 g MgSO <sub>4</sub>	186006812
4 g MgSO <sub>4</sub> , 1 g NaCl, 1 g trisodium citrate dehydrate, 0.5 g disodium hydrogencitrate sesquihydrate, 50/pk	186006813

#### DisQuE Cleanup Tubes (Tube 2)

AOAC Method		
Description	Tube Size	P/N
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> and 50 mg PSA, 100/pk	2 mL	186004572
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> , 50 mg PSA and 50 mg C <sub>18</sub> , 100/pk	2 mL	186004830
DisQuE Tube containing: 900 mg MgSO <sub>4</sub> and 300 mg PSA, 50/pk	15 mL	186008077
DisQuE Tube containing: 900 mg MgSO <sub>4</sub> , 300 mg PSA and 300 mg C <sub>18</sub> , 50/pk	15 mL	186008078
DisQuE Tube containing: 1200 mg MgSO <sub>4</sub> and 400 mg PSA, 50/pk	15 mL	186008072
DisQuE Tube containing: 1200 mg MgSO <sub>4</sub> , 400 mg PSA and 400 mg C <sub>18</sub> , 50/pk	15 mL	186008073
DisQuE Tube containing: 1200 mg MgSO <sub>4</sub> , 400 mg PSA, 400 mg C <sub>18</sub> , and 400 mg GCB, 50/pk	15 mL	186008074

#### DisQuE Cleanup Tubes (Tube 2)

CEN Method		
Description	Tube Size	P/N
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> and 25 mg PSA, 100/pk	2 mL	186004831
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> , 25 mg PSA, and 25 mg C <sub>18</sub> , 100/pk	2 mL	186004832
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> , 25 mg PSA, and 2.5 mg GCB, 100/pk	2 mL	186008076
DisQuE Tube containing: 900 mg MgSO <sub>4</sub> , 150 mg PSA, 50/pk	15 mL	186004833
DisQuE Tube containing: 900 mg MgSO <sub>4</sub> , 150 mg PSA, and 150 mg C <sub>18</sub> , 50/pk	15 mL	186004834

#### DisQuE Cleanup Tubes (Tube 2)

Specialty Cleanup Tubes		
Description	Tube Size	P/N
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> and 50 mg C <sub>18</sub> , 100/pk	2 mL	186008075
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> , 25 mg PSA, 25 mg C <sub>18</sub> , and 7 mg GCB, 100/pk	2 mL	186008071
DisQuE Tube containing: 900 mg MgSO <sub>4</sub> , 450 mg PSA, 300 mg C <sub>18</sub> , and 50 mg GCB, 50/pk	15 mL	186008079
DisQuE Tube containing: 150 mg MgSO <sub>4</sub> , 50 mg PSA, 30 mg C <sub>18</sub> , and 30 mg Alumina-N, 100/pk	2 mL	186008081
DisQuE Tube containing: 750 mg MgSO <sub>4</sub> , 250 mg PSA, 150 mg C <sub>18</sub> , and 150 mg Alumina-N, 50/pk	15 mL	186008080



#### Bulk Sorbents

Description	P/N
Graphitized Carbon Black, 25 g Bottle	186004835
C <sub>18</sub> , 100 g Bottle	WAT035672