




DisQuE EXTRACTION AND CLEANUP TUBES AND POUCHES

DisQuE Extraction and Cleanup Tubes and Pouches are available separately for customized applications and method development. The salts contained in the 50 mL tubes are also available in a pouch format for greater flexibility. The cleanup tubes are available in a standard 2 mL size as well as a 15 mL size for sample enrichment.

Ordering Information

DisQuE Dispersive Sample Preparation Products

Description	P/N
Individual Extraction Tubes (Tube 1)	
50 mL Empty Tube for QuEChERS Extraction ■ 50/pk	186006814
DisQuE 50 mL Tube/ AOAC-Acetate	
 ■ DisQuE 50 mL tube containing: 1.5 g Sodium Acetate and 6 g MgSO ₄ , 100/pk	186004571
DisQuE 50 mL Tube/ CEN-Citrate	
 ■ DisQuE 50 mL tube containing: 1 g trisodium citrate dihydrate, 0.5 g disodium hydrogencitrate sesquihydrate, 1 g NaCl and 4 g MgSO ₄ , 100/pk	186004837

Description	P/N
Individual Extraction Pouch	
DisQuE Pouch	
■ 1.5 g sodium acetate, 6 g MgSO ₄ , 50/pk	186006812
 ■ 4 g MgSO ₄ , 1 g NaCl, 1 g trisodium citrate dehydrate, 0.5 g disodium hydrogencitrate sesquihydrate, 50/pk	186006813

DisQuE Cleanup Tubes (Tube 2)

AOAC Method		
Description	Tube Size	P/N
■ DisQuE Tube containing: 150 mg MgSO ₄ and 50 mg PSA, 100/pk	2 mL	186004572
■ DisQuE Tube containing: 150 mg MgSO ₄ , 50 mg PSA and 50 mg C ₁₈ , 100/pk	2 mL	186004830
■ DisQuE Tube containing: 900 mg MgSO ₄ and 300 mg PSA, 50/pk	15 mL	186008077
■ DisQuE Tube containing: 900 mg MgSO ₄ , 300 mg PSA and 300 mg C ₁₈ , 50/pk	15 mL	186008078
■ DisQuE Tube containing: 1200 mg MgSO ₄ and 400 mg PSA, 50/pk	15 mL	186008072
■ DisQuE Tube containing: 1200 mg MgSO ₄ , 400 mg PSA and 400 mg C ₁₈ , 50/pk	15 mL	186008073
■ DisQuE Tube containing: 1200 mg MgSO ₄ , 400 mg PSA, 400 mg C ₁₈ , and 400 mg GCB, 50/pk	15 mL	186008074

DisQuE Cleanup Tubes (Tube 2)

CEN Method		
Description	Tube Size	P/N
■ DisQuE Tube containing: 150 mg MgSO ₄ and 25 mg PSA, 100/pk	2 mL	186004831
■ DisQuE Tube containing: 150 mg MgSO ₄ , 25 mg PSA, and 25 mg C ₁₈ , 100/pk	2 mL	186004832
■ DisQuE Tube containing: 150 mg MgSO ₄ , 25 mg PSA, and 2.5 mg GCB, 100/pk	2 mL	186008076
■ DisQuE Tube containing: 900 mg MgSO ₄ , 150 mg PSA, 50/pk	15 mL	186004833
■ DisQuE Tube containing: 900 mg MgSO ₄ , 150 mg PSA, and 150 mg C ₁₈ , 50/pk	15 mL	186004834

DisQuE Cleanup Tubes (Tube 2)

Specialty Cleanup Tubes		
Description	Tube Size	P/N
■ DisQuE Tube containing: 150 mg MgSO ₄ and 50 mg C ₁₈ , 100/pk	2 mL	186008075
■ DisQuE Tube containing: 150 mg MgSO ₄ , 25 mg PSA, 25 mg C ₁₈ , and 7 mg GCB, 100/pk	2 mL	186008071
■ DisQuE Tube containing: 900 mg MgSO ₄ , 450 mg PSA, 300 mg C ₁₈ , and 50 mg GCB, 50/pk	15 mL	186008079
■ DisQuE Tube containing: 150 mg MgSO ₄ , 50 mg PSA, 30 mg C ₁₈ , and 30 mg Alumina-N, 100/pk	2 mL	186008081
■ DisQuE Tube containing: 750 mg MgSO ₄ , 250 mg PSA, 150 mg C ₁₈ , and 150 mg Alumina-N, 50/pk	15 mL	186008080



Bulk Sorbents

Description	P/N
Graphitized Carbon Black, 25 g Bottle	186004835
C ₁₈ , 100 g Bottle	WAT035672