









Complete Listing of 12 × 32 mm Vials and Accessories

Screw Cap Vials	Clear	Amber	Max Recovery	300 µL PP	750 µL PP	10 mm Cap Clear	Total Recovery	Amber Max
12 × 32 mm								
Vial Number	17	18	19	20	21	22	23	24
Compatible Systems								
Alliance 2690/2695	▪	▪	—	▪	▪	▪	▪	—
Alliance 2790/2795	▪	▪	▪	▪	▪	—	—	▪
ACQUITY	▪	▪	▪	▪	▪	—	▪	▪
TruView LCMS Certified Combination Packs								
Vial, Cap, and Silicone/PTFE Septum	186005660CV	186005667CV	186005668CV	—	—	—	186005669CV	186005664CV
Vial, Cap, and Pre-slit Silicone/PTFE Septum	186005666CV	186005661CV	186005662CV	—	—	—	186005663CV	186005670CV
LCMS Certified Combination Packs								
Vial, Cap, and Silicone/PTFE Septum	600000751CV	600000752CV	600000749CV	—	—	—	600000750CV	600000754CV
Vial, Cap, and Pre-slit Silicone/PTFE Septum	600000668CV	600000669CV	600000670CV	—	—	—	600000671CV	600000755CV
LC/GC Certified Combination Packs								
Bonded Silicone/PTFE Septum	186000272C	186000846C	186000326C	186002640*	186005220*	WAT270946C	186000384C	186003885C
Combination Deactivated*	186000272DV	186000846DV	186000326DV	—	—	WAT270946DV	186000384DV	—
Bonded Pre-slit Silicone/PTFE Septum	186000307C	186000847C	186000327C	186002639*	186005221*	—	186000385C	186003886C
Combination Deactivated*	186000307DV	186000847DV	186000327DV	—	—	—	186000385DV	—
Combination with PE Septumless Cap	186004132C	186004133C	186004168C	186004112*	—	—	186004167C	—
LC/GC Certified Combination Pack with Cap and PTFE Septum	186007193C	186007194C	186007195C	—	—	—	186007197C	186007196C
Certified Combination Pack with Cap and LB Silicone/PTFE Septum	186007199C	186007200C	186007201C	—	—	—	186007203C	186007202C
Vials Only								
Vials Only	186000273	186000848	186002802	186002626	186005219	WAT063300	186002805	—
Deactivated Vials Only	186000273DV	186000848DV	—	—	—	WAT063300DV	—	—
Injectable Volumes Alliance 2690/2695								
Max	1100 µL	1100 µL	—	280 µL	400 µL	1100 µL	950 µL	—
Residual	750 µL	750 µL	—	20 µL	300 µL	750 µL	9 µL	—
Injectable Volumes Alliance 2790/2795								
Max	1700 µL	1700 µL	1500 µL	290 µL	530 µL	1700 µL	—	1500 µL
Residual	170 µL	170 µL	22 µL	10 µL	170 µL	170 µL	—	22 µL

*Not certified.