

Complete Listing of 12 × 32 mm Vials and Accessories

	Clear	Amber	Max Recovery	300 µL PP	750 µL PP	10 mm Cap Clear	Total Recovery	Amber Max
12 × 32 mm								
<b>Vial Number</b>	17	18	19	20	21	22	23	24
<b>Compatible Systems</b>								
Alliance 2690/2695	▪	▪	—	▪	▪	▪	▪	—
Alliance 2790/2795	▪	▪	▪	▪	▪	—	—	▪
ACQUITY	▪	▪	▪	▪	▪	—	▪	▪
<b>Inserts</b>								
300 µL with Poly Spring	WAT094170(DV)	WAT094170(DV)	—	—	—	WAT094170(DV)	—	—
Max Volume Injection/ Max Residual Volume	230 µL/20 µL	230 µL/20 µL	—	—	—	230 µL/20 µL	—	—
150 µL with Poly Spring	WAT094171(DV)	WAT094171(DV)	—	—	—	WAT094171(DV)	—	—
Max Volume Injection/ Max Residual Volume	144 µL/6 µL	144 µL/6 µL	—	—	—	144 µL/6 µL	—	—
<b>Black Screw Cap for TruView Vials</b>								
PTFE/Silicone Septum	186005826	186005826	186005826	—	—	—	186005826	186005826
Pre-slit PTFE/Silicone Septum	186005827	186005827	186005827	—	—	—	186005827	186005827
<b>Light Blue Screw Cap for LCMS Certified Vials</b>								
PTFE/Silicone Septum	186005828	186005828	186005828	—	—	—	186005828	186005828
Pre-slit PTFE/Silicone Septum	186005829	186005829	186005829	—	—	—	186005829	186005829
<b>Screw Cap and Septum-Silicone/PTFE</b>								
PE Septumless Cap	186004169	186004169	186004169	186004169	186004169	—	186004169	186004169
Blue LectraBond	186000274	186000274	186000274	186000274	186000274	—	186000274	186000274
Red LectraBond	186002129	186002129	186002129	186002129	186002129	—	186002129	186002129
Green LectraBond	186002130	186002130	186002130	186002130	186002130	—	186002130	186002130
White LectraBond	186002456	186002456	186002456	186002456	186002456	—	186002456	186002456
Black Cap with PTFE Septum, 100/pk	186007198	186007198	186007198	186007198	186007198	—	186007198	186007198
<b>Screw Cap and Pre-slit Septum-Silicone/PTFE</b>								
Blue LectraBond	186000305	186000305	186000305	186000305	186000305	—	186000305	186000305
Red LectraBond	186002128	186002128	186002128	186002128	186002128	—	186002128	186002128
Green LectraBond	186002127	186002127	186002127	186002127	186002127	—	186002127	186002127
White LectraBond	186002457	186002457	186002457	186002457	186002457	—	186002457	186002457
<b>For Dissolution System</b>								
Pre-assembled Vial, Cap, and Pre-slit Septum	186000989(DV)	186003455	—	—	—	—	—	—
<b>Storage Cap</b>								
Black Solid 9 mm Cap with Silicone/ PTFE Liner for Sample Storage	186007187	186007187	186007187	186007187	186007187	—	186007187	186007187

When (DV) appears beside a number, a deactivated version of the part can be ordered by adding a DV to the right of the part number.