


SAMPLE PREPARATION WITH PROTEINWORKS FOR LARGE MOLECULE QUANTIFICATION

ProteinWorks™ Sample Preparation Kits enable laboratories to routinely achieve standardized, reproducible, and highly sensitive protein quantification via the surrogate peptide approach.

- Easy sample preparation for a wide variety of proteins including biomarkers and antibodies
- Pre-measured, lot-traceable reagents
- Standardized protocol that has built in flexibility for both 3 and 5 step digestion (reduction & alkylation)
- Reproducible and highly effective protein digests with the ProteinWorks Express Digest Kits
- Achieve your best sensitivity by combining the digest kits with the ProteinWorks μ Elution SPE Cleanup Kit

Take the Variability out of your results



Ordering Information

Description	P/N
ProteinWorks eXpress Digest Kit (for the digestion of immunopurified plasma and serum samples)	176003689
ProteinWorks eXpress Digest Start-Up Kit (eXpress Digest Kit with μ Elution SPE Cleanup Kit and Intact mAb Check Standard)	176003696
ProteinWorks eXpress Direct Digest Kit (for the direct digestion of whole, non-immunopurified plasma and serum samples)	176003688
ProteinWorks eXpress Direct Digest Start-Up Kit (eXpress Direct Digest Kit with μ Elution SPE Cleanup Kit and Intact mAb Check Standard)	176003695
ProteinWorks μ Elution SPE Cleanup Kit (for post-digestion cleanup and sample concentration at the peptide level)	186008304
Intact mAb Check Standard	186006552
Modular Heat Block for 1 mL Tubes	186007985
ACQUITY UPLC Peptide BEH C ₁₈ , 300Å, 1.7 μ m, 2.1 \times 150 mm	186003687
96-well Sample Collection Plate, 700 μ L Round Well	186005837
Polypropylene Cap Mat, Round Well	186002483
96-well Extraction Plate Vacuum Manifold	186001831
Positive Pressure 96-Processor	186006961

Visit www.waters.com/proteinworks to learn more.