


## Solvent Filtration Apparatus

The 300 mL capacity 47 mm Glass Filter Funnel and 1 L capacity 47 mm Glass Funnel/Support Assembly are ideal for vacuum filtration of liquids and degassing of HPLC solvent and mobile phases. The 100% borosilicate glass construction assures resistance to even the most aggressive solvents.

### Ordering Information

#### Solvent Filtration Apparatus

Description	P/N
Solvent Filtration Apparatus 110 V, 60 Hz	176002986
Solvent Filtration Apparatus 220 V, 50 Hz	176002987
All Glass Filter Holder 47 mm, complete	WAT200543
Funnel, 300 mL	WAT200545
Glass Base, tabulated cap	WAT200546
Ground Joint Flask	WAT200547
Swinney Holder	WAT200566
Vacuum Pump 110 V, 60 Hz	725000417
Vacuum Pump 220 V, 50 Hz	725000604



#### Solvent Filtration Membranes

Description	Diameter	Pore Size	Qty.	P/N
PVDF Filter	47 mm	0.45 µm	100/pk	WAT200530
Nylon Filter	47 mm	0.45 µm	100/pk	WAT200532
PTFE Filter	47 mm	0.45 µm	100/pk	WAT200534
	13 mm	0.45 µm	100/pk	WAT200536
GH Polypro Filter	47 mm	0.45 µm	100/pk	WAT200537
Supor (PES) Filter	47 mm	0.45 µm	100/pk	WAT200538
	13 mm	0.45 µm	100/pk	WAT200540
PVDF Filter	47 mm	0.2 µm	100/pk	WAT200531
Nylon Filter	47 mm	0.2 µm	100/pk	WAT200533
PTFE Filter	47 mm	0.2 µm	100/pk	WAT200535
GHP	47 mm	0.2 mm	100/pk	186003524
Supor (PES) Filter	47 mm	0.2 µm	100/pk	WAT200539