









Snap and Crimp Cap Vials	Clear	Amber	Max Recovery	300 µL PP	750 µL PP	Clear Glass Crimp	Amber Crimp	Total Recovery
12 × 32 mm								
Vial Number	25	26	27	28	29	30	31	32
<b>Compatible Systems</b>								
Alliance 2690/2695	▪	▪	—	▪	▪	▪	▪	▪
Alliance 2790/2795	▪	▪	▪	▪	▪	▪	▪	—
ACQUITY	▪	▪	▪	▪	▪	▪	▪	▪
<b>Combination Packs</b>								
Vial, Cap, and Silicone/wPTFE Septum	—	—	—	186002642	186005223	—	—	186000234(DV)
Vial, Cap, and Pre-slit Silicone/PTFE Septum	—	—	—	186002641	186005222	—	—	—
<b>Vials</b>								
Vials Only	WAT094219	WAT094220	186000984	186002628	186005224	WAT094222	WAT094223	186000302
Deactivated Vials Only	WAT094219DV	WAT094220DV	186000984DV	—	—	WAT094222DV	WAT094223DV	186000302DV
<b>Injectable Volumes Alliance 2690/2695</b>								
Max	1100 µL	1100 µL	—	280 µL	400 µL	1100 µL	1100 µL	950 µL
Residual	750 µL	750 µL	—	20 µL	300 µL	750 µL	750 µL	9 µL
<b>Injectable Volumes Alliance 2790/2795</b>								
Max	1700 µL	1700 µL	1500 µL	290 µL	530 µL	1700 µL	1700 µL	—
Residual	170 µL	170 µL	22 µL	10 µL	170 µL	170 µL	170 µL	—
<b>Inserts</b>								
300 µL with Poly Spring	WAT094170(DV)	WAT094170(DV)	—	—	—	WAT094170(DV)	WAT094170(DV)	—
Max Volume Injection/Max Residual Volume	230 µL/20 µL	230 µL/20 µL	—	—	—	230 µL/20 µL	230 µL/20 µL	—
150 µL with Poly Spring	WAT094171(DV)	WAT094171(DV)	—	—	—	WAT094171(DV)	WAT094171(DV)	—
Max Volume Injection/Max Residual Volume	144 µL/6 µL	144 µL/6 µL	—	—	—	144 µL/6 µL	144 µL/6 µL	—
<b>Snap Cap and Septum-Silicone/PTFE</b>								
Blue	186000303	186000303	186000303	186000303	186000303	—	—	186000303
Black	186002649	186002649	186002649	186002649	186002649	—	—	186002649
Red	186002650	186002650	186002650	186002650	186002650	—	—	186002650
<b>Snap Cap and Pre-slit Septum-Silicone/PTFE</b>								
Blue	186000304	186000304	186000304	186000304	186000304	—	—	186000304
Black	186002648	186002648	186002648	186002648	186002648	—	—	186002648
Red	186002647	186002647	186002647	186002647	186002647	—	—	186002647
<b>Snap Cap and PTFE Septum</b>								
Blue	186000328	186000328	186000328	186000328	186000328	—	—	186000328
Black	186002645	186002645	186002645	186002645	186002645	—	—	186002645
Red	186002646	186002646	186002646	186002646	186002646	—	—	186002646
<b>Crimp Cap</b>								
Crimp Cap Silicone/PTFE Septum	—	—	—	—	—	PSL404219	PSL404219	—
Crimp Cap PTFE/Silicone/PTFE Septum	—	—	—	—	—	PSL404231	PSL404231	—
Crimp Cap with Silicone/PTFE Septa	—	—	—	—	—	186006967	186006967	—
Crimper	—	—	—	—	—	PSL904301	PSL904301	—

All items come in quantities of 100 unless otherwise noted.

When (DV) appears beside the part number, a deactivated version of this product can be ordered by adding DV to the right of the part number.