

## Questionnaire for packed GC columns

The following information is necessary to execute you inquiry:

(Please use a separate sheet for each column type)

- |  |  |
|--|--|
| 1. Quantity of the desired column type | pieces   |
| 2. Material of column:                 | Stainless Steel<br>INNO-Steel (silicia-coated stainless steel)<br>other: |

### 3. Dimensions

(Essential for glass columns!

Winding diameter is 160 mm, if other diameter is required please indicate:

mm


The inside diameter is 2 mm as standard and the outer diameter is 1/8"; otherwise please specify!

ID:            and AD:


Length of column:            m

Please mark column geometry type

A



B



7. Packing material	% Coating with: on support: Particle size:		mesh
8. Specifics	Shimadzu type fittings (1/8") (not in general compatible with every shimadzu instrument) Shimadzu: always 5x3 mm Enlargement for injector side end unpacked injector side unpacked detector side		mm mm
9. Fittings	Stainless Steel	Brass	without
10. Conditioning/testing of the columns (only available for stainless steel columns)	without conditioning testing		
11. Comments			